

# NovaGuard® SA Pro Safety System: Summary of Design Verification Testing (DVT)

## I. Introduction

### A. Executive Summary

This technical report summarizes the Design Verification Testing (DVT) of the NovaGuard® SA Pro 1.0 mL Long and 0.5 mL safety systems. The successful results obtained from this testing confirms the safety and efficacy of the NovaGuard® SA Pro 1.0 mL Long and 0.5 mL safety systems design in meeting their design specifications.

### B. Background

The NovaGuard® SA Pro 1.0 mL Long and 0.5 mL safety systems are syringe accessories designed to accompany an ISO 11040-4 compliant 1.0 mL Long and 0.5 mL Pre-filled Syringe (PFS) respectively with a Rigid Needle Shield (RNS) (See Figure 2). Their purpose is to help prevent needle stick injury after the syringe is used. When operated by user action, the NovaGuard® SA Pro safety system elongates to cover the used needle (see Figure 1).

The NovaGuard® SA Pro 1.0 mL long safety system shown in Figure 2 is composed of two plastic injection molded components (Sleeve ①, Syringe Holder ③) and one metal component (Compression Spring ②). The NovaGuard® SA Pro 0.5 mL safety system is identical but shorter in overall length. The NovaGuard® SA Pro safety system integrates a syringe clipping mechanism into the syringe holder (see Figure 2, component ③). The result: syringe insertion becomes a single-step process, enabling existing syringe assembly lines to easily accommodate the NovaGuard® SA Pro safety system with minimal adjustment and change parts.





Figure 1: How it works



## NovaGuard® SA Pro Safety System: Summary of Design Verification Testing (DVT)

---

### II. Experimental Details of Design Verification Testing (DVT)

1. NovaGuard® SA Pro safety system test samples were subject to climate conditioning (as per ASTM D4332-14); and transport simulation testing (as per ASTM D4169-16, ASTM D5276-98, ASTM D999-08 and ASTM D4728-06) at an independent test facility prior to DVT functional testing.
2. One batch was used for T0 (unaged) testing, with the remaining batches placed into temperature and humidity-controlled environments for aging. Batches were stored at accelerated aging conditions and tested at specified time intervals to simulate up to 5 years real-time equivalent; the remaining batches are being stored at room temperature conditions and will be tested each year, real-time, up to five years.
3. The below list describes key tests performed during Design Verification Testing.
  - Transport Simulation
  - Syringe assembly force
  - Activation security and sound
  - Pre-activation disassembly force
  - Drop test, pre- and post-activation
  - Activation and safety lock
  - Tensile disassembly, pre- and post-activation
  - Safety force, post-activation
  - Needle retraction depth, post-activation
  - RNS removal and replacement

### III. Summary

***Unaged batches (T0) for NovaGuard® SA Pro 1.0 mL Long and NovaGuard® SA Pro 0.5 mL safety systems passed each test by meeting the pre-determined acceptance criteria.***

***Accelerated aged batches (up to 5 years real-time equivalent) for NovaGuard® SA Pro 1.0 mL Long and NovaGuard® SA Pro 0.5 mL safety systems passed each test by meeting the pre-determined acceptance criteria.***

***The results obtained from design verification testing of unaged and accelerated aged NovaGuard® SA Pro safety system confirms the adequacy of the NovaGuard® SA Pro safety system in meeting the design input specifications and a shelf life of 5 years.***

***Multi-year real-time aging and testing of NovaGuard® SA Pro 1.0 mL Long and NovaGuard® SA Pro 0.5 mL safety systems up to 5 years is on-going.***

## NovaGuard® SA Pro Safety System: Summary of Design Verification Testing (DVT)

### IV. Acronyms

Acronym	Explanation
ASTM	American Society for Testing and Materials
DVT	Design Verification Testing
ISO	International Standards Organisation
PFS	Pre-Filled Syringe
RNS	Rigid Needle Shield

### V. Reference Documents

Document Number	Title
ISO 11040-4/A1: 2020	Pre-Filled Syringes – Part 4: Glass barrels for injectables
EN ISO 23908: 2013	Sharps Injury Protection – Requirements and Test Methods
EN ISO 11608-1: 2015	Needle-based Injection systems for medical use
ASTM D4169-16	Standard Practice for Performance Testing of Shipping Containers and Systems
ASTM D4332-14	Practice for Conditioning Containers, Packages, or Packaging Components for Testing
ASTM D999-08 (2015)	Test Methods for Vibration Testing of Shipping Containers
ASTM D4728-06 (2012)	Test Method for Random Vibration Testing of Shipping Containers
ASTM D5276-98 (2009)	Test Method for Drop Test of Loaded Containers by Free Fall

This document is for informational purposes only. West's products and solutions are sold on the basis that it is the customer's responsibility to evaluate and test the West product or solution to determine its compatibility with other materials and fitness for any end use. WEST MAKES NO WARRANTIES, WHETHER EXPRESS, IMPLIED OR STATUTORY, INCLUDING, WITHOUT LIMITATION, THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR PARTICULAR PURPOSE, RELATING TO THE INFORMATION IN THIS DOCUMENT.

This technical report dated 28 March 2022 is the second release version of this report and supersedes the previous version.

West and the diamond logo and NovaGuard® are registered trademarks of West Pharmaceutical Services, Inc. in the United States and other jurisdictions.

©2022 by West Pharmaceutical Services, Inc.

All rights reserved. This material is protected by copyright. No part of it may be reproduced, stored in a retrieval system, or transmitted in any form or by any means, electronic, mechanical, photocopying or otherwise, without written permission of West Pharmaceutical Services, Inc.