

InsoCap®



The High-Performance Choice for Blow/Fill/Seal IV Bottles

InsoCap closures for IV bottles provide a latex-free solution for Blow/Fill/Seal operations. West's unique two-stage process molds the plastic shell and thermoplastic elastomer septum in a continuous operation. The injection area is sealed with foil to protect against contamination and tampering.

InsoCap benefits:

- Delivered ready-to-use
- Closures fit existing filling lines and Eurohead bottles
- Protective foil, available in gold and silver, can be decorated for branding
- No leakage between the septum and the plastic cap
- Easily recyclable

InsoCap attributes:

- Shell is available in transparent or white polypropylene and polyethylene
- Septum is moulded from West's proprietary TPE formulation
- Closures comply with ISO 15759, specific to plastic caps for Blow/Fill/Seal processing
- Orientation feature positions the cap for ease of penetration

InsoCap provides superior performance for:

- Seal integrity between the closure and the IV bottle
- Retaining the inserted IV spike
- Creating a tight seal between the TPE septum and the IV spike or needle
- Resealability
- IV spike or needle removal force
- Ease of removing the foil seal

To learn more about InsoCap closures, contact your West account manager or a technical customer service representative at TCS.Europe-AsiaPacific@westpharma.com.



westpharma.com

Denmark
Horsens
+45 (7561) 6000

England
St. Austell
+44 (1726) 6 35 63

France
Saint Germain en Laye
+33 (1) 39 21 54 00

Germany
Eschweiler
+49 (24 03) 7960

Italy
Milan
+39 02-6703079

Serbia
Beograd
+38 11 132 85 430

Spain
Madrid
+34 (91) 656-1415

West warrants that its products will conform to all agreed-upon specifications or those specifications published by West. No one is authorised to make any warranty or representation other than as set forth in this standard warranty, and no customer or user may rely on any other warranty or representation. All recommendations and suggestions contained herein regarding end uses and methods of use are based upon laboratory studies, test results and experience available to West. However, West cannot control the conditions or circumstances under which the company's products may be used, and accordingly, all users have sole responsibility to evaluate and test West's products for compatibility and fitness for use. WEST MAKES NO OTHER WARRANTIES, EXPRESS OR IMPLIED, AND DISCLAIMS ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR PARTICULAR PURPOSE. West's responsibility under its warranty is limited to the cost of replacing West's products sold to the original purchaser. A West account representative can provide more information and answer specific questions.
InsoCap® is a registered trademark of West Pharmaceutical Services, Inc.