



Leading the way with integrated containment and delivery solutions

SMARTDOSE® 10 Injector

PRECISE. RELIABLE.
ADAPTIVE.

Expanded Platform

- Leveraging the commercial success of SmartDose 3.5 injector
- Designed to handle large volumes and variety of delivery times through adaptable, pre-programmable technology
- Proven engineering minimizes risks & industrialization minimizes investment and time to market

Wearable Injector

- Making treatment more patient-friendly with home administration
- Ability to deliver high volume and high viscosity drug products
- Integrated connectivity enables home-based patient monitoring from the clinic

Patient-Centric

- User-experience, easy to use, intuitive design
- Visual, tactile and audible feedback to boost user confidence
- Onboarding and training solutions

Integrated Solutions for Fill Finish

- West can enable fill-finish and simplify your path to market
- Small scale, in-house lab filling services
- GMP filling with established CMOs
- West program management and analytical services to assist new line build & retrofits

 **PRECISE**

 **RELIABLE**

 **ADAPTIVE**



Technical Specifications	
Target Tissue	Subcutaneous
Max Container Volume	Up to 10mL
Sterilization	Ready to Use sterile
Container Type	Custom Daikyo Crystal Zenith® cartridge with septum/cap & Flurotec® film-coated piston
Container Material	Daikyo Crystal Zenith cyclic olefin polymer (COP) USP
Tub & Nest Format	10mL: 5x6 configuration for standard filling equipment
Needle Safety System	Complies with ISO standard
Driving Mechanism	Electromechanical
System Control	Pre-programmable
Cartridge Loading	User loaded
Wireless Connectivity	Yes
Temperature Sensing	Optional capability

For more information visit www.westpharma.com/SmartDose

North America 1-800-345-9800 option 8 | **South America** +55 11 4055 6061 | **Europe** +49 (0) 2403-7960 | **Asia Pacific** +65 6862 3400

West Pharmaceutical Services, Inc. 530 Herman O. West Drive | Exton, PA 19341