




## VIAL2BAG ADVANCED™

*Universal, Needle-Free Admixture  
Device for Immediate Use  
IV Administration*

The Vial2Bag Advanced™ 20mm admixture device enables reconstitution and transfer of a drug between a vial and an IV bag prior to administration to the patient. The innovative needle-free and easy-to-use design is ideal for immediate use and can optimize pharmacy clean room time, while assisting in standardizing admixture procedures to maximize efficiency of point of care nursing.

### New Design

- Dual channel design provides dedicated fluid pathways into and out of the IV bag
- Enhanced vial spike design improves the connection between vial and Vial2Bag Advanced™ 20mm adapter and the transfer of medication
- Needle-free device meets OSHA and NAPPSI sharps prevention guidance



### Optimized Patient Care

- Provides a simple solution for just-in-time IV medication preparation, improving timeliness at point of care
- Easy to implement immediate use solution

### Effective Waste Reduction

- Reduces medication and IV fluid waste
- Universal compatibility with all manufacturers' 50, 100, and 250mL IV bags
- Suitable for use with liquid and powder medications with standard 20mm closure vials



**Vial2Bag Advanced™ 20mm admixture device can be an effective solution when meeting the compounding provisions associated with USP <797>.\***

\*Note: The use of the Vial2Bag Advanced™ 20mm admixture device in association with any provision of compounding should not be interpreted as modifying, extending or superseding a drug manufacturer's labeling recommendations for storage and expiration dating. Refer to drug manufacturer's recommendations and USP compounding guidelines for shelf life and sterility information



**Vial2Bag Advanced™  
20mm Admixture Device**  
Item # 36098150

**Quantity**  
50/Box  
150/Carton



**Distributed by: Progressive Medical, Inc.** 997 Horan Dr. • Fenton, MO 63026  
**Ordering Information:** Phone: 314-961-5786 • Fax: 314-961-4535  
cs@progressivemedinc.com



## Innovation in Reconstitution and Transfer Systems

Scan to learn more about our products and services

For more information visit [www.westpharma.com](http://www.westpharma.com)

**North America** 1-800-345-9800 option 8 | **South America** +55 11 4055 6061 | **Europe** +49 (0) 2403-7960 | **Asia Pacific** +65 6862 3400

**West Pharmaceutical Services, Inc.** 530 Herman O. West Drive | Exton, PA 19341