

## Leading the way with integrated containment and delivery solutions



# SMARTDOSE®

GEN. III – 3.5mL

PRECISE. RELIABLE.  
PRELOADED.

### Expanded Platform

- Pre-loaded SmartDose Generation (Gen.) III system
- Uses proprietary prefillable cartridge as primary container
- Final assembly does not require cleanroom environment

### Wearable Injector

- Reduced user steps with pre-loaded cartridge
- Subcutaneous injection for self-administration
- Pre-programmable injection time and rate
- Ability to deliver high volume and high viscosity drug products

### Patient-Centric

- Simple, easy-to-use system with fewer user steps
- Visual, tactile and audible feedback boosts patient confidence
- Onboarding and training solutions



For more information visit [www.westpharma.com](http://www.westpharma.com)

-  **PRECISE**
-  **RELIABLE**
-  **PRELOADED**



### Technical Specifications

Target Tissue	Subcutaneous
Max Container Volume	Up to 3.7mL
Sterilization	Device sterilization not required Ready to Use (RU) sterile container system
Container Type	Custom cartridge with cannula/needle shield and FluroTec® piston
Container Material	Daikyo Crystal Zenith® cyclic olefin polymer (COP)
Tub & Nest Format	8x4 configuration for standard filling equipment
Needle Safety System	Automatic needle insertion and retraction
Driving Mechanism	Electromechanical
System Control	Pre-programmable
Cartridge Loading	Pre-loaded
Wireless Connectivity	Bluetooth
Temperature Sensing	Optional capability

West Pharmaceutical Services, Inc. | 530 Herman O. West Drive, Exton, PA 19341 | [www.smartdose.com](http://www.smartdose.com) | North America +1 800-345-9800  
Europe +49 2403 7960  
Asia Pacific +64 6862 3400

West seeks partners for its SmartDose drug delivery platform. This platform is intended to be used as an integrated system with drug filling and final assembly completed by the pharmaceutical/biotechnology company.

West and the diamond logo, By your side for a healthier world and FluroTec® are registered trademarks or trademarks of West Pharmaceutical Services, Inc., in the United States and other jurisdictions. SmartDose® is a registered trademark of Medimop Medical Projects Ltd., a subsidiary of West Pharmaceutical Services, Inc. Daikyo Crystal Zenith® is a registered trademark of Daikyo Seiko, Ltd. Daikyo Crystal Zenith® and FluroTec® technologies are licensed from Daikyo Seiko, Ltd.