

# Leading the way with integrated containment and delivery solutions



## SMARTDOSE®

PRECISE. RELIABLE.  
READY TO GO.

### FDA Approved

- The first combination product that incorporates the SmartDose platform technology was recently approved by the US Food and Drug Administration (FDA)
- Thousands of doses have been administered using the SmartDose platform
- Proven engineering, manufacturing and regulatory expertise to support your needs

### Wearable Injector

- Subcutaneous self-administration
- Ability to deliver high volume and high viscosity drug products

### Patient-Centric

- User-centered design
- Connectivity to a variety of software platforms
- Able to link with adherence solutions
- Onboarding and training solutions available

### Flexible Technology

- Address a variety of delivery times through adaptable, pre-programmable technology
- Maximize patient comfort through pre-programmable delivery times

In the US, the SmartDose platform is considered an on-body infusor with prefilled cartridge. This drug delivery system incorporates human factors and usability to deliver a truly patient-centric approach.

For more information about the SmartDose platform and how it can be used to enhance delivery of your drug product, contact your account manager.

-  **PRECISE**
-  **RELIABLE**
-  **READY TO GO**

## TECHNICAL SPECIFICATIONS

Target Tissue	Subcutaneous
Max Container Fill Volume	Up to 3.7mL
Sterilization	Ready to use Sterile
Container Type	Daikyo Crystal Zenith®, a proprietary cyclic olefin polymer (COP), cartridge is capped
Tub & Nest Format	8x8 nested
Needle Safety System	Passive needle safety system
Driving Mechanism	Electromechanical
System Control	Micro-processor and firmware controlled
Cartridge Loading	User loaded



West Pharmaceutical Services, Inc. | 530 Herman O. West Drive, Exton, PA 19341 | [www.smartdose.com](http://www.smartdose.com) | North America +1 800-345-9800  
 Europe +49 2403 7960  
 Asia Pacific +65 6860 5879

West seeks partners for its SmartDose platform. This platform is intended to be used as an integrated system with drug filling and final assembly completed by the pharmaceutical/biotechnology company. West and the diamond logo, By your side for a healthier world and FluroTec® are registered trademarks or trademarks of West Pharmaceutical Services, Inc., in the United States and other jurisdictions. SmartDose® is a registered trademark of Medimop Medical Projects Ltd., a subsidiary of West Pharmaceutical Services, Inc. Daikyo Crystal Zenith® is a registered trademark of Daikyo Seiko, Ltd. Daikyo Crystal Zenith® and FluroTec® technologies are licensed from Daikyo Seiko, Ltd.