



By your side
for a healthier world™

Leading the way with integrated containment and delivery solutions



SMARTDOSE® 10 Injector

**PRECISE. RELIABLE.
ADAPTIVE.**

➤ Expanded Platform

- Leveraging the commercial success of the SmartDose 3.5 injector
- Designed to handle large volumes and various delivery times through adaptable, pre-programmable technology
- Proven engineering minimizes risks & industrialization minimizes investment and time to market

➤ Wearable Injector

- Making treatment more patient-friendly with home administration
- Ability to deliver high volume and high viscosity drug products
- Integrated connectivity enables home-based patient monitoring from the clinic

➤ Patient-Centric

- User-experience, easy to use, intuitive design
- Visual, tactile and audible feedback to boost user confidence
- Onboarding and training solutions

➤ Integrated Solutions for Fill Finish

- West can enable fill-finish and simplify your path to market
- Small scale, in-house lab filling services
- GMP filling with established CMOs
- West program management and analytical services to assist new line build & retrofits

 **PRECISE**

 **RELIABLE**

 **ADAPTIVE**



Technical Specifications

Target Tissue	Subcutaneous
Max Container Volume	Up to 10mL
Sterilization	Ready to Use sterile
Container Type	Custom Daikyo Crystal Zenith® cartridge with septum/cap & Flurotec® film-coated piston
Container Material	Daikyo Crystal Zenith cyclic olefin polymer (COP) USP
Tub & Nest Format	10mL: 5x6 configuration for standard filling equipment
Needle Safety System	Complies with ISO standard
Driving Mechanism	Electromechanical
System Control	Pre-programmable
Cartridge Loading	User loaded
Wireless Connectivity	Yes
Temperature Sensing	Optional capability

For more information visit www.westpharma.com/SmartDose

North America 1-800-345-9800 option 8 | South America +55 11 4055 6061 | Europe +49 (0) 2403-7960 | Asia Pacific +65 6862 3400

West Pharmaceutical Services, Inc. 530 Herman O. West Drive | Exton, PA 19341

West and the diamond logo and By your side for a healthier world are registered trademarks or trademarks of West Pharmaceutical Services, Inc., in the United States and other jurisdictions. SmartDose and logo are trademarks of West Pharma. Services IL, Ltd., a subsidiary of West Pharmaceutical Services, Inc., in the United States and other jurisdictions. Crystal Zenith and Flurotec are registered trademark of Daikyo Seiko, Ltd. Daikyo Crystal Zenith technology is licensed from Daikyo Seiko, Ltd.