


- Notes:
1. Drawing may not be to scale.
 2. Max. radius of R0.15[0.006] allowed on sharp corners.
 3. Values denoted by [X.XXX] are in inches.


UNCONTROLLED COPY

For internal use only



SECTION A-A

N/A	N/A	N/A	N/A	N/A	N/A
N/A	N/A	N/A	N/A	N/A	N/A
0	N/A	N/A	SAP 200048626	N/A	N/A
Rev.	DR.	Rev. Date	Description	App.	App. Date
Revision Table					

Surface Area: 1.12cm² [.174in²]			
Tolerance : ISO 3302-1 Class M 3			
Units: Millimeters			
Approvals	Date	Title: 1 mL long plug NovaPure plunger	
DWN: S.Chajar	21 Sep 2023		
CHKD: T.Eck	27 Sep 2023		
Approved H.Giang	27 Sep 2023	Article: 1 mL long plug	
Reference Drawings		Material: 4023/50 Grey	
X-6476-P1 Rev F	N/A	Drawing No.	Revision
N/A	N/A	WP-465	0
First Angle Projection		Scale: 10:1	Sheet Size: A3
		The information on this drawing is proprietary. This information shall not be disclosed, reproduced in any form or used for any purpose, without prior written consent of West Pharmaceutical Services.	