

4

3

2

1

PRIMARY BAG –
STERILIZABLEBAG – Heat Sealed

SECONDARY BAG –
STERILIZABLEBAG – Heat Sealed

TERTIARY CLEAR
CR100 PE BAG – Heat Sealed

PLASTIC CARTON
15 1/4" X 12 1/4" X 2 1/4" O.D.
(1) DOUBLE STERILIZABLEBAG/CARTON

CARTON LABEL

LABEL ID


LABEL ID

LABEL ID

USE BY DATE

STERILE
CLAIM



THIS DRAWING IS FOR INFORMATION ONLY AND IS NOT INTENDED FOR CONTROLLED SPECIFICATION RELEASE.					
TOLERANCES UNLESS NOTED		APPROVALS	DATE		
DECIMAL	METRIC	OWN		MATERIAL:	
.XX ± .01	.X ± .3	A. Barrilli	23 JUN 16	DRAWING AND REV. NO:	
.XXX ± .003	.XX ± .08	M. Habash	23 JUN 16	1794-144	
.XXXX ± .0015	.XXX ± .038	M. Mariani	23 JUN 16	SCALE: 1:1	
FRACTIONAL INCH: ± 1/64		TOOL ENGR.		WEIGHT: LB.	
ANGULAR: ± 1°		PREC. ENGR.		SHEET 1 OF 1	
SURFACE FINISH: ✓		APPD.		SIZE: C	
THE INFORMATION ON THIS DRAWING IS PROPRIETARY. THIS INFORMATION SHALL NOT BE DISCLOSED, REPRODUCED IN ANY FORM OR USED FOR ANY PURPOSE, WITHOUT PRIOR WRITTEN CONSENT OF WEST PHARMACEUTICAL SERVICES.				SAP Number	
DO NOT SCALE DRAWING				TOOL CODE:	