



- Drawing may be not to scale
- Radii without dimension 0.15 - 0.2 [0.006 - 0.008]

- Nominal surface area = 12,31 cm² [1,908 inch²]



Rubber Formulation: 4023/50 grey

Revision Nb. 4

Unlimited dimensions ISO 3302, Part 1 Class M 3

Replaced by: 1346 from 25/Apr/08

Scale: 5:1	drawn: Hei	17/Jun/20	approved:
------------	------------	-----------	-----------

17/Jun/20

17/Jun/20

FluroTec® - Lyo Stopper

Article Nb.:
1346

THE INFORMATION ON THIS DRAWING IS PROPRIETARY.
THIS INFORMATION SHALL NOT BE DISCLOSED, REPRODUCED IN ANY FORM OR USED FOR ANY PURPOSE,
WITHOUT PRIOR WRITTEN CONSENT OF WEST PHARMACEUTICAL SERVICES.