

Packaging Options That Fit Your Needs

Our commitment to the entire commercial life cycle of a drug product means we give our component packaging just as much time and attention as we give the components themselves. We've done the work to ensure our components are easily integrated into your manufacturing operations.

From concept to market, we are there by your side, helping to meet the many packaging requirements you need for your components.




SMALL VOLUME

- Ideal for testing or sampling needs
- Available products ship within 72 hours
- Westar® washed, gamma or steam sterilized bag formats



BULK BAGS

- Ideal for traditional cleanroom filling lines
- Applicable for use with restricted access barrier system (RABS) filling lines
- Westar® washed, gamma or steam sterilized bag formats



PORTED BAGS

- Designed for use with alpha-beta rapid transfer port (RTP) systems
- Interlocking system designed to open only when both port doors are properly docked and integrated
- Ideal for restricted access barrier system (RABS) and isolator filling lines
- Westar® washed, gamma or steam sterilized bag formats

West Pharmaceutical Services, Inc. | 530 Herman O. West Drive, Exton, PA 19341 | www.westpharma.com

West and the diamond logo, By your side for a healthier world™ and Westar® are registered trademarks of West Pharmaceutical Services, Inc., in the United States and other jurisdictions.

For complete contact information please visit www.westpharma.com.

Copyright © 2015 West Pharmaceutical Services, Inc.

#9120 0615